



fkpA Antibody, Biotin conjugated

Product Code	CSB-PA517042HD01AXW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O08437
Immunogen	Recombinant Aeromonas hydrophila FKBP-type peptidyl-prolyl cis-trans isomerase FkpA protein (21-268AA)
Raised In	Rabbit
Species Reactivity	Aeromonas hydrophila
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	FKBP-type peptidyl-prolyl cis-trans isomerase FkpA (PPIase) (EC 5.2.1.8) (Rotamase), fkpA
Immunogen Species	Aeromonas hydrophila
Research Area	Others
Target Names	fkpA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.