



DAG1 Antibody, Biotin conjugated

Product Code	CSB-PA613486ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14118
Immunogen	Peptide sequence from Human Dystroglycan protein (726-744AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dystroglycan (Dystrophin-associated glycoprotein 1) [Cleaved into: Alpha-dystroglycan (Alpha-DG); Beta-dystroglycan (Beta-DG)], DAG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DAG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.