



GRIK2 Antibody, Biotin conjugated

Product Code	CSB-PA618751LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13002
Immunogen	Recombinant Human Glutamate receptor ionotropic, kainate 2 protein (203-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor ionotropic, kainate 2 (GluK2) (Excitatory amino acid receptor 4) (EAA4) (Glutamate receptor 6) (GluR-6) (GluR6), GRIK2, GLUR6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.