



STK4 Antibody, HRP conjugated

Product Code	CSB-PA621644NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13043
Immunogen	Peptide sequence from Human Serine/threonine-protein kinase 4 protein (6-24AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase 4 (EC 2.7.11.1) (Mammalian STE20-like protein kinase 1) (MST-1) (STE20-like kinase MST1) (Serine/threonine-protein kinase Krs-2) [Cleaved into: Serine/threonine-protein kinase 4 37kDa subunit (MST1/N); Serine/threonine-protein kinase 4 18kDa subunit (MST1/C)], STK4, KRS2 MST1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.