



LAMA4 Antibody, HRP conjugated

Product Code	CSB-PA621948LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16363
Immunogen	Recombinant Human Laminin subunit alpha-4 protein (1593-1823AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laminin subunit alpha-4 (Laminin-14 subunit alpha) (Laminin-8 subunit alpha) (Laminin-9 subunit alpha), LAMA4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LAMA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.