



RASA2 Antibody, HRP conjugated

Product Code	CSB-PA622992LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15283
Immunogen	Recombinant Human Ras GTPase-activating protein 2 protein (609-850AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras GTPase-activating protein 2 (GTPase-activating protein 1m) (GAP1m), RASA2, GAP1M RASGAP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.