



RASA2 Antibody, Biotin conjugated

Product Code	CSB-PA622992LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15283
Immunogen	Recombinant Human Ras GTPase-activating protein 2 protein (609-850AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras GTPase-activating protein 2 (GTPase-activating protein 1m) (GAP1m), RASA2, GAP1M RASGAP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASA2