



FAM110C Antibody, HRP conjugated

Product Code	CSB-PA637287LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1W6H9
Immunogen	Recombinant Human Protein FAM110C protein (50-321AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM110C, FAM110C
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM110C