



FAM110C Antibody, Biotin conjugated

Product Code	CSB-PA637287LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1W6H9
Immunogen	Recombinant Human Protein FAM110C protein (50-321AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM110C, FAM110C
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM110C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.