



NKAPL Antibody, Biotin conjugated

Product Code	CSB-PA702797ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5M9Q1
Immunogen	Peptide sequence from Human NKAP-like protein (114-123AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NKAP-like protein, NKAPL, C6orf194
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NKAPL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.