



OTUD3 Antibody, Biotin conjugated

Product Code	CSB-PA719399LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T2D3
Immunogen	Recombinant Human OTU domain-containing protein 3 protein (108-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	OTU domain-containing protein 3 (EC 3.4.19.12), OTUD3, KIAA0459
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OTUD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.