



PLA2G4F Antibody, FITC conjugated

Product Code	CSB-PA721207LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68DD2
Immunogen	Recombinant Human Cytosolic phospholipase A2 zeta protein (74-299AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic phospholipase A2 zeta (cPLA2-zeta) (EC 3.1.1.4) (Phospholipase A2 group IVF), PLA2G4F
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLA2G4F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.