



PPM1L Antibody, HRP conjugated

Product Code	CSB-PA722509LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SGD2
Immunogen	Recombinant Human Protein phosphatase 1L protein (128-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1L (EC 3.1.3.16) (Protein phosphatase 1-like) (Protein phosphatase 2C isoform epsilon) (PP2C-epsilon), PPM1L, PP2CE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPM1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.