



DAB2IP Antibody, HRP conjugated

Product Code	CSB-PA730138OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VWQ8
Immunogen	Peptide sequence from Human Disabled homolog 2-interacting protein (85-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disabled homolog 2-interacting protein (DAB2 interaction protein) (DAB2-interacting protein) (ASK-interacting protein 1) (AIP-1) (DOC-2/DAB-2 interactive protein), DAB2IP, AF9Q34 AIP1 KIAA1743
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	DAB2IP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.