



CAPRIN2 Antibody, HRP conjugated

Product Code	CSB-PA740238OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IMN6
Immunogen	Peptide sequence from Human Caprin-2 protein (83-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caprin-2 (C1q domain-containing protein 1) (Cytoplasmic activation/proliferation-associated protein 2) (Gastric cancer multidrug resistance-associated protein) (Protein EEG-1) (RNA granule protein 140), CAPRIN2, C1QDC1 EEG1 KIAA1873 RNG140
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CAPRIN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.