



CAPRIN2 Antibody, FITC conjugated

Product Code	CSB-PA740238OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IMN6
Immunogen	Peptide sequence from Human Caprin-2 protein (83-100AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Caprin-2 (C1q domain-containing protein 1) (Cytoplasmic activation/proliferation-associated protein 2) (Gastric cancer multidrug resistance-associated protein) (Protein EEG-1) (RNA granule protein 140), CAPRIN2, C1QDC1 EEG1 KIAA1873 RNG140
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CAPRIN2