



ENHO Antibody, HRP conjugated

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA744265OB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6UWT2 |
| Immunogen | Peptide sequence from Human Adropin protein (37-54AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Adropin (Energy homeostasis-associated protein), ENHO, C9orf165 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | ENHO |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |