



TNFAIP8L2 Antibody, HRP conjugated

Product Code	CSB-PA754193LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P589
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 8-like protein 2 protein (1-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor alpha-induced protein 8-like protein 2 (TIPE2) (TNF alpha-induced protein 8-like protein 2) (TNFAIP8-like protein 2) (Inflammation factor protein 20), TNFAIP8L2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TNFAIP8L2