



TNFAIP8L2 Antibody, Biotin conjugated

Product Code	CSB-PA754193LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P589
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 8-like protein 2 protein (1-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor alpha-induced protein 8-like protein 2 (TIPE2) (TNF alpha-induced protein 8-like protein 2) (TNFAIP8-like protein 2) (Inflammation factor protein 20), TNFAIP8L2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TNFAIP8L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.