



DNAJC10 Antibody, HRP conjugated

Product Code	CSB-PA809006LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXB1
Immunogen	Recombinant Human DnaJ homolog subfamily C member 10 protein (521-793AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily C member 10 (EC 1.8.4.-) (Endoplasmic reticulum DNA J domain-containing protein 5) (ER-resident protein ERdj5) (ERdj5) (Macrothioredoxin) (MTHr), DNAJC10, ERDJ5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJC10