



DNAJC10 Antibody, Biotin conjugated

Product Code	CSB-PA809006LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXB1
Immunogen	Recombinant Human DnaJ homolog subfamily C member 10 protein (521-793AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily C member 10 (EC 1.8.4.-) (Endoplasmic reticulum DNA J domain-containing protein 5) (ER-resident protein ERdj5) (ERdj5) (Macrothioredoxin) (MTHr), DNAJC10, ERDJ5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJC10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.