



TNS4 Antibody, FITC conjugated

Product Code	CSB-PA810319LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZW8
Immunogen	Recombinant Human Tensin-4 protein (20-231AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tensin-4 (C-terminal tensin-like protein), TNS4, CTEN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNS4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.