





SUPV3L1 Antibody

Product Code	CSB-PA814232LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYB8
Immunogen	Recombinant Human ATP-dependent RNA helicase SUPV3L1, mitochondrial protein (499-786AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent RNA helicase SUPV3L1, mitochondrial (EC 3.6.4.13) (Suppressor of var1 3-like protein 1) (SUV3-like protein 1), SUPV3L1, SUV3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SUPV3L1

Image



IHC image of CSB-PA814232LA01HU diluted at 1:300 and staining in paraffin-embedded human prostate tissue performed on a Leica BondTM system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.


Western Blot

Positive WB detected in: Hela whole cell lysate,
K562 whole cell lysate

All lanes: SUPV3L1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 88 kDa

Observed band size: 88 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.