



SUPV3L1 Antibody, Biotin conjugated

Product Code	CSB-PA814232LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYB8
Immunogen	Recombinant Human ATP-dependent RNA helicase SUPV3L1, mitochondrial protein (499-786AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent RNA helicase SUPV3L1, mitochondrial (EC 3.6.4.13) (Suppressor of var1 3-like protein 1) (SUV3-like protein 1), SUPV3L1, SUV3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SUPV3L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.