



TMEM237 Antibody, HRP conjugated

Product Code	CSB-PA822287OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96Q45
Immunogen	Peptide sequence from Human Transmembrane protein 237 protein (59-78AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 237 (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 4 protein), TMEM237, ALS2CR4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM237
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.