



AJUBA Antibody, HRP conjugated

Product Code	CSB-PA836229NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96IF1
Immunogen	Peptide sequence from Human LIM domain-containing protein ajuba protein (61-77AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM domain-containing protein ajuba, AJUBA, JUB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AJUBA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.