



DTX3 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA839833LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8N9I9 |
| Immunogen | Recombinant Human Probable E3 ubiquitin-protein ligase DTX3 protein (90-347AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Probable E3 ubiquitin-protein ligase DTX3 (EC 2.3.2.27) (Protein deltex-3) (Deltex3) (RING finger protein 154) (RING-type E3 ubiquitin transferase DTX3), DTX3, RNF154 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | DTX3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |