



DTX3 Antibody, Biotin conjugated

Product Code	CSB-PA839833LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9I9
Immunogen	Recombinant Human Probable E3 ubiquitin-protein ligase DTX3 protein (90-347AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable E3 ubiquitin-protein ligase DTX3 (EC 2.3.2.27) (Protein deltex-3) (Deltex3) (RING finger protein 154) (RING-type E3 ubiquitin transferase DTX3), DTX3, RNF154
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DTX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.