



WFIKKN1 Antibody, HRP conjugated

Product Code	CSB-PA842729OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NZ8
Immunogen	Peptide sequence from Human WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1 protein (187-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1 (Growth and differentiation factor-associated serum protein 2) (GASP-2) (hGASP-2) (WAP, follistatin, immunoglobulin, Kunitz and NTR domain-containing protein), WFIKKN1, C16orf12 GASP2 WFIKKN
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WFIKKN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.