



WFIKKN1 Antibody, FITC conjugated

Product Code	CSB-PA842729OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NZ8
Immunogen	Peptide sequence from Human WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1 protein (187-207AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1 (Growth and differentiation factor-associated serum protein 2) (GASP-2) (hGASP-2) (WAP, follistatin, immunoglobulin, Kunitz and NTR domain-containing protein), WFIKKN1, C16orf12 GASP2 WFIKKN
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WFIKKN1