



LEO1 Antibody, HRP conjugated

Product Code	CSB-PA845150NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVC0
Immunogen	Peptide sequence from Human RNA polymerase-associated protein LEO1 protein (2-24AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA polymerase-associated protein LEO1 (Replicative senescence down-regulated leo1-like protein), LEO1, RDL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	LEO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.