



LEO1 Antibody, Biotin conjugated

Product Code	CSB-PA845150ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVC0
Immunogen	Peptide sequence from Human RNA polymerase-associated protein LEO1 protein (2-24AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA polymerase-associated protein LEO1 (Replicative senescence down-regulated leo1-like protein), LEO1, RDL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	LEO1