



RANBP9 Antibody, HRP conjugated

Product Code	CSB-PA853501OB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96S59
Immunogen	Peptide sequence from Human Ran-binding protein 9 protein (158-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ran-binding protein 9 (RanBP9) (BPM-L) (BPM90) (Ran-binding protein M) (RanBPM) (RanBP7), RANBP9, RANBPM
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RANBP9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.