



RANBP9 Antibody, Biotin conjugated

Product Code	CSB-PA853501OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96S59
Immunogen	Peptide sequence from Human Ran-binding protein 9 protein (158-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ran-binding protein 9 (RanBP9) (BPM-L) (BPM90) (Ran-binding protein M) (RanBPM) (RanBP7), RANBP9, RANBPM
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RANBP9