



ELMO1 Antibody, FITC conjugated

Product Code	CSB-PA856404HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92556
Immunogen	Recombinant Human Engulfment and cell motility protein 1 protein (520-727AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Engulfment and cell motility protein 1 (Protein ced-12 homolog), ELMO1, KIAA0281
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ELMO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.