



BBC3 Antibody, HRP conjugated

Product Code	CSB-PA861160NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXH1
Immunogen	Peptide sequence from Human Bcl-2-binding component 3 protein (5-23AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bcl-2-binding component 3 (JFY-1) (p53 up-regulated modulator of apoptosis), BBC3, PUMA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BBC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.