



BBC3 Antibody, FITC conjugated

Product Code	CSB-PA861160NC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXH1
Immunogen	Peptide sequence from Human Bcl-2-binding component 3 protein (5-23AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bcl-2-binding component 3 (JFY-1) (p53 up-regulated modulator of apoptosis), BBC3, PUMA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BBC3