



WDR83 Antibody, Biotin conjugated

Product Code	CSB-PA863110LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRX9
Immunogen	Recombinant Human WD repeat domain-containing protein 83 protein (64-315AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat domain-containing protein 83 (Mitogen-activated protein kinase organizer 1) (MAPK organizer 1), WDR83, MORG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	WDR83
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.