



MIOS Antibody, HRP conjugated

Product Code	CSB-PA865167LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NXC5
Immunogen	Recombinant Human GATOR complex protein MIOS protein (1-237AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	GATOR complex protein MIOS (Missing oocyte meiosis regulator homolog), MIOS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MIOS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.