




DDX47 Antibody

Product Code	CSB-PA867123LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0S4
Immunogen	Recombinant Human Probable ATP-dependent RNA helicase DDX47 protein (225-460AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable ATP-dependent RNA helicase DDX47 (EC 3.6.4.13) (DEAD box protein 47), DDX47
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DDX47

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, HL-60 whole cell lysate, K562 whole cell lysate, 293T whole cell lysate

All lanes: DDX47 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 51, 46 kDa

Observed band size: 51 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.