



DDX47 Antibody, HRP conjugated

Product Code	CSB-PA867123LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0S4
Immunogen	Recombinant Human Probable ATP-dependent RNA helicase DDX47 protein (225-460AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable ATP-dependent RNA helicase DDX47 (EC 3.6.4.13) (DEAD box protein 47), DDX47
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DDX47
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.