



DAXX Antibody, HRP conjugated

Product Code	CSB-PA871395LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UER7
Immunogen	Recombinant Human Death domain-associated protein 6 protein (77-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Death domain-associated protein 6 (Daxx) (hDaxx) (ETS1-associated protein 1) (EAP1) (Fas death domain-associated protein), DAXX, BING2 DAP6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DAXX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.