



AIFM2 Antibody, Biotin conjugated

Product Code	CSB-PA874794LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRQ8
Immunogen	Recombinant Human Apoptosis-inducing factor 2 protein (110-373AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptosis-inducing factor 2 (EC 1.-.-) (Apoptosis-inducing factor homologous mitochondrion-associated inducer of death) (Apoptosis-inducing factor-like mitochondrion-associated inducer of death) (p53-responsive gene 3 protein), AIFM2, AMID PRG3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	AIFM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.