



# AIFM2 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA874794LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BRQ8
<b>Immunogen</b>	Recombinant Human Apoptosis-inducing factor 2 protein (110-373AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Apoptosis-inducing factor 2 (EC 1.-.-) (Apoptosis-inducing factor homologous mitochondrion-associated inducer of death) (Apoptosis-inducing factor-like mitochondrion-associated inducer of death) (p53-responsive gene 3 protein), AIFM2, AMID PRG3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	AIFM2