



TJP2 Antibody, Biotin conjugated

Product Code	CSB-PA880071ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UDY2
Immunogen	Peptide sequence from Human Tight junction protein ZO-2 protein (1165-1183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tight junction protein ZO-2 (Tight junction protein 2) (Zona occludens protein 2) (Zonula occludens protein 2), TJP2, X104 ZO2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TJP2