



SH2D3A Antibody, Biotin conjugated

Product Code	CSB-PA880088OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRG2
Immunogen	Peptide sequence from Human SH2 domain-containing protein 3A protein (211-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH2 domain-containing protein 3A (Novel SH2-containing protein 1), SH2D3A, NSP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SH2D3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.