



TCF7L1 Antibody, Biotin conjugated

Product Code	CSB-PA884626OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCS4
Immunogen	Peptide sequence from Human Transcription factor 7-like 1 protein (80-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor 7-like 1 (HMG box transcription factor 3) (TCF-3), TCF7L1, TCF3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TCF7L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.