



WWC3 Antibody, Biotin conjugated

Product Code	CSB-PA890929ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULE0
Immunogen	Peptide sequence from Human Protein WWC3 protein (65-85AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein WWC3, WWC3, KIAA1280
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WWC3