



WWC3 Antibody, Biotin conjugated

Product Code	CSB-PA890929ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULE0
Immunogen	Peptide sequence from Human Protein WWC3 protein (65-85AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein WWC3, WWC3, KIAA1280
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WWC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.