



PYCARD Antibody, Biotin conjugated

Product Code	CSB-PA890936LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULZ3
Immunogen	Recombinant Human Apoptosis-associated speck-like protein containing a CARD protein (1-195AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptosis-associated speck-like protein containing a CARD (hASC) (Caspase recruitment domain-containing protein 5) (PYD and CARD domain-containing protein) (Target of methylation-induced silencing 1), PYCARD, ASC CARD5 TMS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PYCARD